

Šūnu kultūru reaģenti

APRAKSTS

ScienCell piedāvā pilnu komplektu augstas kvalitātes šūnu kultūru reaģentu

Produkti:

- [0.2% Gelatin Solution, 100 ml](#)
- [100X Non-Essential Amino Acids, 100 ml](#)
- [3D Collagen I Matrix Kit](#)
- [Antibiotic/Antimycotic Solution, 50 ml](#)
- [Antimycotic Solution, 50 ml](#)
- [Bovine Hypothalamus Extract \(BHE-1\)](#)
- [Bovine Hypothalamus Extract \(BHE-5\)](#)
- [Bovine Pituitary Extract, 100 mg](#)
- [Bovine Pituitary Extract, 25 mg](#)
- [Bovine Plasma Fibronectin](#)
- [Bovine Plasma Vitronectin](#)
- [Cell Culture Grade Water, 500 ml](#)
- [Cell Freezing Medium, 50 ml](#)
- [Cell Freezing Medium-serum free, 50 ml](#)
- [Dulbecco's Phosphate-Buffered Saline, 500 ml](#)
- [Enzyme-free Cell Dissociation Solution, 100 ml](#)
- [Fetal Bovine Serum, 10 ml](#)
- [Fetal Bovine Serum, 25 ml](#)
- [Fetal Bovine Serum, 500 ml](#)
- [Glucose Solution, 100ml](#)
- [Hank's Balanced Salt Solution](#)
- [Heparin \(0.2%\), 5 ml](#)
- [Human Cellular Fibronectin](#)
- [Human Pluripotent Stem Cell Cryopreservation Medium \(5 x 10 ml\)](#)
- [Human Pluripotent Stem Cell Cryopreservation Medium, 50 ml](#)
- [Insulin Transferrin Selenium, 10 ml](#)
- [Iron Fortified Bovine Serum](#)
- [L-Glutamine Solution, 100 ml](#)
- [Mycoplasma Elimination Solution \(10 mg/ml; 1 ml\)](#)
- [Penicillin/Streptomycin Solution, 100 ml](#)
- [Penicillin/Streptomycin Solution](#)
- [Poly-L-Lysine, 1 mg/ml](#)
- [Poly-L-Lysine, 10 mg/ml](#)
- [Puromycin \(10 mg/ml; 10 ml\)](#)
- [Puromycin \(10mg/ml; 1 ml\)](#)
- [Soybean Trypsin Inhibitor](#)
- [StemDS® Human Pluripotent Stem Cell Dissociation Solution](#)
- [Trypan Blue Solution \(0.4%\)](#)
- [Trypsin Neutralization Solution, 100 ml](#)
- [Trypsin Neutralization Solution, 200 ml](#)
- [Trypsin/EDTA \(0.05%; 200ml\)](#)
- [Trypsin/EDTA \(0.25%; 200ml\)](#)
- [Trypsin/EDTA, 0.05%](#)
- [Trypsin/EDTA, 0.25%](#)

