

## P-30/100, platform

### DESCRIPTION

Platform with clamps for flasks.

### SPECIFICATIONS

|                  |              |
|------------------|--------------|
| Working area     | 360 × 400 mm |
| Flask volume     | 100-150 ml   |
| Number of clamps | 30           |



### CAT. NUMBER

BS-010135-BK

Platform P-30/100

### ACCESSORIES



**ES-20/60**  
Without platform  
Orbital Shaker-Incubator

Orbital Shaker-Incubator **ES-20/60** for biotechnological and pharmaceutical laboratories is a professional category equipment designed for cultivation of microorganisms and eukaryotic cells including animal, plant and insect cells. It is also ...

[read more](#)



**PSU-20i**  
Without platform  
Multi-functional Orbital Shaker

Shaker **PSU-20i** provides three motion types: orbital, reciprocal and vibrating, which can be performed separately, pairwise and sequentially in repeated cycles. Shaker is designed for applications both in small specialized ...

[read more](#)