


# HSC-1000, clamp

## DESCRIPTION

Tight fit clamp 1000 ml - Ø130 mm for UP-168



## CAT. NUMBER

BS-010167-IK

HSC-1000, clamp 1000ml

## ACCESSORIES



**UP-168**  
BS-010135-JK  
platform

Universal platform with clamps can accommodate flasks or bottles of different volume sizes. Clamps are not included with platform and need to be ordered separately.



**ES-20/80**  
Software included, without  
platform  
Orbital Shaker-Incubator

ES-20/80 shaker-incubator for biotechnological and pharmaceutical laboratories is a professional category equipment. The typical applications include - microbial and cell culture cultivation, protein expression, solubility studies, general mixing, as well ...

[read more](#)



**ES-20/80C**  
Software included, without  
platform  
Orbital Shaker-Incubator with  
cooling

ES-20/80C shaker-incubator for biotechnological and pharmaceutical laboratories is a professional category equipment. The typical applications include - microbial and cell culture cultivation, protein expression, solubility studies, general mixing, as well ...

[read more](#)