

## P-2/24, platform

### DESCRIPTION

For 24 tubes 1.5-2 ml.

### SPECIFICATIONS

Capacity	24
Microtes tubes volume	1.5-2 ml



### CAT. NUMBER

BS-010216-AK

Platform P-2/24



#### MPS-1

With built-in universal platform  
High-Speed Multi Plate Shaker

High-Speed Multi Plate Shaker **MPS-1** can be used in virtually any application by providing adjustable mixing of reagents in microtest plates, PCR plates, deepwell plates and test tubes (shaking tubes ...

[read more](#)