

## AP-96, adapter

### DESCRIPTION

2 adapters for 96-well semi-skirted and non-skirted PCR plates - made of Ertacetal® C. Autoclavable



### CAT. NUMBER

BS-010219-DK

Adapter set (2 pcs) - AP-96



**LMC-3000**  
Without rotor  
Laboratory Centrifuge

**LMC-3000** is a modern low-speed bench-top centrifuge designed for operation with microtest plates and centrifuge tubes up to 50 ml. This device is widely used in biomedical profile laboratories.

Features: ...

[read more](#)



**LMC-12**  
Without rotor  
Laboratory Centrifuge

**LMC-3000** is a modern low-speed bench-top centrifuge designed for operation with microtest plates and centrifuge tubes up to 50 ml. This device is widely used in biomedical profile laboratories.

Features: ...

[read more](#)



**CVP-2**  
Including adapters AP-96  
Centrifuge-Vortex for PCR plates

After many years of Combined Centrifuge-Vortex concept success, we are proud to introduce the long awaited Centrifuge/Vortex for PCR plates, **CVP-2**, to the sample preparation market.

The Spin-Mix-Spin technology ...

[read more](#)



**LMC-4200R**  
Without rotor  
Refrigerated Laboratory Centrifuge

Laboratory bench-top centrifuge with refrigeration **LMC-4200R** provides temperature control of biomaterial during centrifugation. Temperature control of the so-called "cold-shelf" is a gold standard for enzymologists and cell biologists because it ...

[read more](#)



**LMC-56**  
Without rotor  
Laboratory Centrifuge

**LMC-56** is a modern benchtop, low-speed laboratory centrifuge designed for convenient sedimentation, centrifugation and collection of necessary samples. It provides operation with tubes, blood collection systems, gel cards, microtest and ...

[read more](#)



**R-2**  
BS-010208-AK  
Rotor

Swing-out rotor for standard 96-well microtitre plates, skirted PCR plates, deepwell plates up to 45mm.



**Ri-2**  
BS-010221-GK  
Rotor

Swing-out rotor for standard 96-well microtitre plates, skirted PCR plates, deepwell plates up to 45mm.