

## TR-16/19, rack

### DESCRIPTION

Test tube rack for  $\varnothing 16$  to  $\varnothing 19$  mm tubes, capacity 16 tubes. WxDxH: 155x90x112 mm

### SPECIFICATIONS

Capacity

16

Tube size

$\varnothing 16-19$  mm



### CAT. NUMBER

BS-010406-FK

Test tube rack for  $\varnothing 16$  to  $\varnothing 19$  mm tubes, capacity 16 tubes



**WB-4MS**  
With base BP-1 and lid  
Stirred water bath

Water bath-thermostat **WB-4MS** is designed for chemical, pharmaceutical, medical and biological laboratory research.

**WB-4MS** provides increased temperature stabilization (up to 0.1°C) due to built-in magnetic stirrer (regulated speed 250–1,000 rpm).

...

[read more](#)