

# HAS-1, Hand operator set

## DESCRIPTION

HAS-1 Hand operator set

- 1ch. adapter
  - with ejector for 200 mkl tips
  - with 40 mm stainless steel pin
  - with 125 mm stainless steel pin
- 8ch. adapter
  - with ejector for 200 mkl tips
  - with 35 mm stainless steel pins
- Stand SH-6



## CAT. NUMBER

BS-040118-PK

HAS-1, Hand operator set



**FTI-2i**  
BS-040120-A02  
Aspirator with Trap Flask

Aspirator with trap flask **FTI-2i** is designed for aspiration or removal of alcohol, buffer and liquid from reaction vessels (e.g. during DNA/RNA purification or other macromolecule reprecipitation techniques).

The device ...

[read more](#)