


External pre-filter set (carbon/1µm) with manometer

DESCRIPTION



CAT. NUMBER

BS-070104-KK

External pre-filter set (carbon/1µm) with manometer

ACCESSORIES



Labaqua Trace
BS-070105-A02
ultrapure water system

Labaqua ultrapure systems are multi-purpose water purification systems. The Labaqua systems produce ultrapure and pure water directly from tap water.

Any configuration of a Labaqua ultrapure system produces both ultrapure ...

[read more](#)



Labaqua HPLC
BS-070104-A02
ultrapure water system

Labaqua ultrapure systems are multi-purpose water purification systems. The Labaqua systems produce ultrapure and pure water directly from tap water.

Labaqua HPLC produces water with very low organic carbon (TOC) ...

[read more](#)



Labaqua Bio
30l storage tank
ultrapure water system

Labaqua ultrapure systems are multi-purpose water purification systems. The Labaqua systems produce ultrapure and pure water directly from tap water.

Labaqua Bio system produces water with very low organic and ...

[read more](#)