


UV bulb 185 nm

DESCRIPTION



CAT. NUMBER

BS-070104-DK

UV bulb 185 nm

ACCESSORIES



Labagua Trace
BS-070105-A02
ultrapure water system

Labagua ultrapure systems are multi-purpose water purification systems. The Labagua systems produce ultrapure and pure water directly from tap water.

Any configuration of a Labagua ultrapure system produces both ultrapure ...

[read more](#)



Labagua HPLC
BS-070104-A02
ultrapure water system

Labagua ultrapure systems are multi-purpose water purification systems. The Labagua systems produce ultrapure and pure water directly from tap water.

Labagua HPLC produces water with very low organic carbon (TOC) ...

[read more](#)



Labagua Bio
30l storage tank
ultrapure water system

Labagua ultrapure systems are multi-purpose water purification systems. The Labagua systems produce ultrapure and pure water directly from tap water.

Labagua Bio system produces water with very low organic and ...

[read more](#)