


# Extended tubing

## DESCRIPTION

Extended tubing 2 m long, with fittings and MA-U adapter



## CAT. NUMBER

**BS-040120-DK**

Extended tubing 2 m long, with fittings and MA-U adapter



**FTA-2i**  
BS-040120-A02  
Aspirator with Trap Flask

Aspirator with trap flask **FTA-2i** is designed for aspiration or removal of alcohol, buffer and liquid from reaction vessels (e.g. during DNA/RNA purification or other macromolecule reprecipitation techniques).

The device ...

[read more](#)