


Calibration kit, d12mm

DESCRIPTION

McFarland Turbidity Standards, Ø12mm

SPECIFICATIONS

Diameter	12 mm
----------	-------



CAT. NUMBER

21255	ø12mm
-------	-------

ACCESSORIES



DEN-1

BS-050102-AAF

Densitometer (suspension
turbidity detector)

Densitometer is designed for measurement of cell suspension's turbidity in the range of 0.0–6.0 McFarland units ($0 - 180 \times 10^7$ cells/ml).

Densitometer provides the opportunity to measure solution ...

[read more](#)



DEN-1B

BS-050104-AAF

Densitometer (suspension
turbidity detector)

Densitometer is designed for measurement of cell suspension's turbidity in the range of 0.0–6.0 McFarland units ($0 - 180 \times 10^7$ cells/ml).

Densitometers provide the opportunity to measure solution turbidity in ...

[read more](#)



DEN-600

BS-050109-AAA

Photometer

DEN-600 is a compact, portable, rechargeable battery powered photometer. DEN-600 comprises of 600 nm wavelength optical system, which enables to apply - 1) OD600 method that estimated total number of ...

[read more](#)



DEN-600m

IVD certified

Photometer IVD

DEN-600m (IVD version) is a compact, portable, rechargeable battery powered photometer. DEN-600m comprises of 600 nm wavelength optical system, which enables to apply - 1) OD600 method that estimated total ...

[read more](#)