

# RMT-24, Rotor

## DESCRIPTION

Fixed-angle rotor for 24 - 1.5 or 2 ml microtubes

## SPECIFICATIONS

|            |          |
|------------|----------|
| Capacity   | 24       |
| Max. speed | 6000 rpm |
| Max. RCF   | 3750 × g |



## CAT. NUMBER

BS-010221-BK

Rotor RMT-24

## ACCESSORIES



### LMC-4200R

Without rotor  
 Refrigerated Laboratory  
 Centrifuge

Laboratory bench-top centrifuge with refrigeration **LMC-4200R** provides temperature control of biomaterial during centrifugation. Temperature control of the so-called "cold-shelf" is a gold standard for enzymologists and cell biologists because it ...

[read more](#)



### LMC-56

Without rotor  
 Laboratory Centrifuge

**LMC-56** is a modern benchtop, low-speed laboratory centrifuge designed for convenient sedimentation, centrifugation and collection of necessary samples. It provides operation with tubes, blood collection systems, gel cards, microtest and ...

[read more](#)



### RR-U

BS-010208-UK  
 rack for rotors

Rack for rotors