

## BI-25-6, adapter set

### DESCRIPTION

For rotors R-6, R-6P, Ri-6, Ri-6P optional insert set (6 pcs) - BI-25-6 for 25 ml conical tubes



### CAT. NUMBER

**BS-010221-VK**

For rotors R-6, R-6P, Ri-6, Ri-6P optional insert  
set (6 pcs) - BI-25-6 for 25 ml conical tubes



**LMC-3000**  
Without rotor  
Laboratory Centrifuge

**LMC-3000** is a modern low-speed bench-top centrifuge designed for operation with microtest plates and centrifuge tubes up to 50 ml. This device is widely used in biomedical profile laboratories.

Features: ...

[read more](#)



**LMC-12**  
Without rotor  
Laboratory Centrifuge

**LMC-3000** is a modern low-speed bench-top centrifuge designed for operation with microtest plates and centrifuge tubes up to 50 ml. This device is widely used in biomedical profile laboratories.

Features: ...

[read more](#)



**LMC-4200R**  
Without rotor  
Refrigerated Laboratory Centrifuge

Laboratory bench-top centrifuge with refrigeration **LMC-4200R** provides temperature control of biomaterial during centrifugation. Temperature control of the so-called "cold-shelf" is a gold standard for enzymologists and cell biologists because it ...

[read more](#)



**LMC-56**  
Without rotor  
Laboratory Centrifuge

**LMC-56** is a modern benchtop, low-speed laboratory centrifuge designed for convenient sedimentation, centrifugation and collection of necessary samples. It provides operation with tubes, blood collection systems, gel cards, microtest and ...

[read more](#)



**R-6**  
BS-010208-DK  
rotor

Swing-out rotor for conical bottom centrifuge tube. Aluminium buckets.



**R-6P**  
BS-010208-XK  
Rotor

Swing-out rotor for conical bottom centrifuge tube. Plastic buckets (POM Kocetal).



**Ri-6**  
BS-010221-HK  
rotor

Swing-out rotor for conical bottom centrifuge tube. Aluminium buckets.



**Ri-6P**  
BS-010221-IK  
Rotor

Swing-out rotor for conical bottom centrifuge tube. Plastic buckets (POM Kocetal).