


## B-12AB, block

### DESCRIPTION

Block for one deep-well plate Abgene™ 96/1200 µl - Cat.No. AB0564, AB0787



### CAT. NUMBER

BS-010159-SK

Block for one deep-well plate Abgene™ 96/1200  
µl



### TS-DW

Without thermoblock

Thermo-Shaker for deep well plates

TS-DW Thermo-Shaker is designed for shaking and incubating deep well plates.

A multisystem principle, used in the design of the Thermo-Shaker, allows operating it as 3 independent devices: Incubator, Plate

...

[read more](#)