


B-2PO, block

DESCRIPTION

Block for one deep-well plate PlateOne® 96/2000 µl - Cat.No. S1896-2000



CAT. NUMBER

BS-010159-NK

Block for one deep-well plate PlateOne®
96/2000 µl



TS-DW

Without thermoblock

Thermo-Shaker for deep well plates

TS-DW Thermo-Shaker is designed for shaking and incubating deep well plates.

A multisystem principle, used in the design of the Thermo-Shaker, allows operating it as 3 independent devices: Incubator, Plate

...

[read more](#)