

Remote grade 2 water dispenser

DESCRIPTION

Universal remote dispenser set with 3 m supply hose and water distribution module



CAT. NUMBER

BS-070104-JK

Remote dispenser

ACCESSORIES



Labaqua Trace
BS-070105-A02
ultrapure water system

Labaqua ultrapure systems are multi-purpose water purification systems. The Labaqua systems produce ultrapure and pure water directly from tap water.

Any configuration of a Labaqua ultrapure system produces both ultrapure ...

[read more](#)



Labaqua HPLC
BS-070104-A02
ultrapure water system

Labaqua ultrapure systems are multi-purpose water purification systems. The Labaqua systems produce ultrapure and pure water directly from tap water.

Labaqua HPLC produces water with very low organic carbon (TOC) ...

[read more](#)



Labaqua Bio
30l storage tank
ultrapure water system

Labaqua ultrapure systems are multi-purpose water purification systems. The Labaqua systems produce ultrapure and pure water directly from tap water.

Labaqua Bio system produces water with very low organic and ...

[read more](#)