

Storage tank 60l

DESCRIPTION

Storage tank with base, tap and multipoint level switch, 60 l



CAT. NUMBER

BS-070102-SK

Storage tank, 60l

ACCESSORIES



Labaqua Trace
BS-070105-A02
ultrapure water system

Labaqua ultrapure systems are multi-purpose water purification systems. The Labaqua systems produce ultrapure and pure water directly from tap water.

Any configuration of a Labaqua ultrapure system produces both ultrapure ...

[read more](#)



Labaqua HPLC
BS-070104-A02
ultrapure water system

Labaqua ultrapure systems are multi-purpose water purification systems. The Labaqua systems produce ultrapure and pure water directly from tap water.

Labaqua HPLC produces water with very low organic carbon (TOC) ...

[read more](#)



Labaqua Bio
30l storage tank
ultrapure water system

Labaqua ultrapure systems are multi-purpose water purification systems. The Labaqua systems produce ultrapure and pure water directly from tap water.

Labaqua Bio system produces water with very low organic and ...

[read more](#)