

T-4L, Table

DESCRIPTION

New modular design of laboratory furniture provides flexibility and ease of use.

SPECIFICATIONS

Maximum load	75 kg
Overall dimensions (W×D×H)	1290 × 600 × 770 mm
Weight	36 kg



CAT. NUMBER

BS-040107-BK

T-4L, Table

ACCESSORIES



UVT-S-AR
BS-040107-AAA
DNA/RNA UV-cleaner box

DNA/RNA UV-cleaner box **UVT-S-AR** is designed for clean operations with DNA samples. UV-cleaner box provide protection against contamination.

Model is a bench-top type, made of metal framework, glass walls, working ...

[read more](#)



LF-1
BS-050101-BK
Drawer unit

New modular design of laboratory furniture provides flexibility and ease of use.