

## P-6/250, platform

### DESCRIPTION

Platform with clamps for flasks.

### SPECIFICATIONS

Working area	250 × 190 mm
Flask volume	250-300 ml
Number of clamps	6



### CAT. NUMBER

BS-010108-DK

Platform P-6/250

### ACCESSORIES



**ES-20**  
 Without platform  
 Orbital Shaker-Incubator

The **ES-20** is a compact bench-top Shaker-Incubator used for mixing of biological liquids as well as for incubation and cultivation of biological liquids according to the operator set program

Built-in ...

[read more](#)



**PSU-10i**  
 Without platform  
 Orbital Shaker

The **PSU-10i** provides regulated orbital motion of the platform and is designed for use both in small specialized biotechnological laboratories and in large multidisciplinary laboratories: a choice of five (5) ...

[read more](#)



**CPS-20**  
 Without platform  
 CO<sub>2</sub> Shaker

The **CPS-20** provides regulated orbital motion of the platform and is designed for use specifically in CO<sub>2</sub> incubators. CPS-20 is specifically designed for use in harsh environments such as ...

[read more](#)