

P-9/500, platform

DESCRIPTION

Platform with clamps for flasks.

SPECIFICATIONS

Working area	360 × 400 mm
Flask volume	500 ml
Number of clamps	9



CAT. NUMBER

BS-010135-AK

Platform P-9/500

ACCESSORIES



ES-20/60
 Without platform
 Orbital Shaker-Incubator

Orbital Shaker-Incubator **ES-20/60** for biotechnological and pharmaceutical laboratories is a professional category equipment designed for cultivation of microorganisms and eukaryotic cells including animal, plant and insect cells. It is also ...

[read more](#)



PSU-20i
 Without platform
 Multi-functional Orbital Shaker

Shaker **PSU-20i** provides three motion types: orbital, reciprocal and vibrating, which can be performed separately, pairwise and sequentially in repeated cycles. Shaker is designed for applications both in small specialized ...

[read more](#)