

TR-21/50, rack

DESCRIPTION

Adjustable angle test tube rack for UP-168

SPECIFICATIONS

| | |
|-----------------|-------|
| Tube volume | 50 ml |
| Number of tubes | 21 |
| Hole diameter | 30 mm |



CAT. NUMBER

BS-010135-KK

TR-21/50, test tube rack, 50ml

ACCESSORIES



UP-168
BS-010135-JK
platform

Universal platform with clamps can accommodate flasks or bottles of different volume sizes. Clamps are not included with platform and need to be ordered separately.



ES-20/60
Without platform
Orbital Shaker-Incubator

Orbital Shaker-Incubator **ES-20/60** for biotechnological and pharmaceutical laboratories is a professional category equipment designed for cultivation of microorganisms and eukaryotic cells including animal, plant and insect cells. It is also ...

[read more](#)



ES-20/80
Software included, without
platform
Orbital Shaker-Incubator

ES-20/80 shaker-incubator for biotechnological and pharmaceutical laboratories is a professional category equipment. The typical applications include - microbial and cell culture cultivation, protein expression, solubility studies, general mixing, as well ...

[read more](#)



ES-20/80C
Software included, without
platform
Orbital Shaker-Incubator with
cooling

ES-20/80C shaker-incubator for biotechnological and pharmaceutical laboratories is a professional category equipment. The typical applications include - microbial and cell culture cultivation, protein expression, solubility studies, general mixing, as well ...

[read more](#)



PSU-20i
Without platform
Multi-functional Orbital Shaker

Shaker **PSU-20i** provides three motion types: orbital, reciprocal and vibrating, which can be performed separately, pairwise and sequentially in repeated cycles. Shaker is designed for applications both in small specialized ...

[read more](#)