

PP-400, platform

DESCRIPTION

Flat platform with non-slip silicone mat can accommodate various low profile containers.

SPECIFICATIONS

Working area

360 × 400 mm



CAT. NUMBER

BS-010135-FK

Platform PP-400 with silicone mat

ACCESSORIES



ES-20/60

Without platform
 Orbital Shaker-Incubator

Orbital Shaker-Incubator **ES-20/60** for biotechnological and pharmaceutical laboratories is a professional category equipment designed for cultivation of microorganisms and eukaryotic cells including animal, plant and insect cells. It is also ...

[read more](#)



ES-20/80

Software included, without platform
 Orbital Shaker-Incubator

ES-20/80 shaker-incubator for biotechnological and pharmaceutical laboratories is a professional category equipment. The typical applications include - microbial and cell culture cultivation, protein expression, solubility studies, general mixing, as well ...

[read more](#)



ES-20/80C

Software included, without platform
 Orbital Shaker-Incubator with cooling

ES-20/80C shaker-incubator for biotechnological and pharmaceutical laboratories is a professional category equipment. The typical applications include - microbial and cell culture cultivation, protein expression, solubility studies, general mixing, as well ...

[read more](#)