

PRS-48, platform

DESCRIPTION

48 tubes 10-16 mm diameter (1.5 ml-15 ml tubes).

SPECIFICATIONS

Platform type	standard
Tube diameter	10-16 mm
Tube volume	1.5-15 ml
Capacity	48



CAT. NUMBER

BS-010118-CK

Platform PRS-48



Multi RS-60
Including PRS-48 platform
Programmable rotator

Programmable Rotator **Multi RS-60** performs several motion types in one module. Microprocessor control allows performing not only Vertical rotation (1) of the platform, but also Reciprocal rotation (2) as well ...

[read more](#)