

# SC-18/02C, block

## DESCRIPTION

20 × 0.2 ml microtubes + 12 × 1.5 ml microtubes

## SPECIFICATIONS

Amount of 0.2 ml tubes	20
Amount of 1.5 ml tubes	12



## CAT. NUMBER

BS-010143-CK	SC-18/02C
--------------	-----------

## ACCESSORIES



TS-100C

Without thermoblock

Thermo-Shaker with cooling for microtubes and PCR plates

Thermo-Shaker TS-100C provides intensive mixing and temperature control of samples in microtest tubes or PCR plate. This model of Thermo-Shaker differs from TS-100 with a possibility of cooling samples down ...

[read more](#)



TS-100C Smart

Software included, without thermoblock

Thermo-Shaker with cooling for microtubes and PCR plates

Thermo-Shaker TS-100C Smart provides intensive mixing and temperature control of samples in microtest tubes or PCR plate. This model of Thermo-Shaker differs from TS-100 with a possibility of cooling samples ...

[read more](#)