

## R-1.5, rotor

### DESCRIPTION

Rotor for 12 x 1.5/2 ml microtest tubes

### SPECIFICATIONS

Amount of 2 ml tubes	12
Amount of 1.5 ml tubes	12



### CAT. NUMBER

BS-010205-AK	Rotor R-1.5
--------------	-------------

### ACCESSORIES



#### FVL-2400N

Including rotors R-1.5, R-0.5/0.2  
Minicentrifuge-Vortex Combi-Spin

Minicentrifuge-vortex Combi-Spin **FVL-2400N** is specially designed for genetic engineering research (for PCR-diagnostics experiments). Units can be used in microbiological, biochemical, clinical laboratories and industrial biotechnological laboratories.

Combi-Spin provides simultaneous mixing...

[read more](#)



#### MSC-3000

Including rotors R-1.5, R-0.5/0.2  
Minicentrifuge-Vortex Multi-Spin

Minicentrifuge-Vortex Multi-Spin **MSC-3000** is product of many years evolution of Spin-Mix-Spin technology that is intended for collecting micro volumes of reagents on the microtube's bottom (first centrifugation spin), following mixing...

[read more](#)



#### MSC-6000

Including rotors R-1.5, R-0.5/0.2  
Minicentrifuge-Vortex Multi-Spin

Minicentrifuge-Vortex Multi-Spin **MSC-6000** is product of many years evolution of Spin-Mix-Spin technology that is intended for collecting micro volumes of reagents on the microtube's bottom (first centrifugation spin), following mixing...

[read more](#)