

A-02, adapter

■ DESCRIPTION

Adapter 12 pieces for microtubes 0.2 ml

■ SPECIFICATIONS

Amount of 0.2 ml tubes

12



■ CAT. NUMBER

BS-010213-BK A-02 adapter - 12 pieces for microtubes 0.2 ml

**Microspin 12**

Including MSL-SC lid, adapters A-02, A-05
High-speed Mini-centrifuge

High-speed mini-centrifuge **Microspin 12** is a compact desktop centrifuge designed for biomedical laboratories.

Microspin 12 is used for extracting RNA/DNA samples, sedimentation of biological components, biochemical and chemical analysis of

...

[read more](#)