

A-05, adapter

DESCRIPTION

Adapter 12 pieces for microtubes 0.5 ml

SPECIFICATIONS

Amount of 0.5 ml tubes

12



CAT. NUMBER

BS-010213-AK A-05 adapter - 12 pieces for microtubes 0.5 ml



Microspin 12

Including MSL-SC lid, adapters A-02, A-05
High-speed Mini-centrifuge

High-speed mini-centrifuge **Microspin 12** is a compact desktop centrifuge designed for biomedical laboratories.

Microspin 12 is used for extracting RNA/DNA samples, sedimentation of biological components, biochemical and chemical analysis of

...

[read more](#)