


USB 2.0 Hub 10 × ports

DESCRIPTION



CAT. NUMBER

BS-010158-BK

USB 2.0 Hub 10 × ports

ACCESSORIES



RTS-1

Including TPP TubeSpin®
 Bioreactor vessels 50ml, 20pcs
 Personal bioreactor

RTS-1 is personal bioreactor which utilize patented Reverse-Spin® technology that applies non-invasive, mechanically driven, low energy consumption, innovative type of agitation where cell suspension is mixed by the singleuse falcon ...

[read more](#)



RTS-1C

Including TPP TubeSpin®
 Bioreactor vessels 50ml, 20pcs
 Personal bioreactor

RTS-1C is personal bioreactor which utilize patented Reverse-Spin® technology that applies non-invasive, mechanically driven, low energy consumption, innovative type of agitation where cell suspension is mixed by the singleuse falcon ...

[read more](#)



RTS-8

Including TPP TubeSpin®
 Bioreactor vessels 50ml, 20pcs
 Personal multi-channel bioreactor with non-invasive real time OD measurement

RTS-8 is a personal bioreactor that utilizes patented Reverse-Spin® technology that applies non-invasive, mechanically driven, low energy consumption, innovative type of agitation where cell suspension is mixed by the single-use ...

[read more](#)



RTS-8 Plus

Including TPP TubeSpin®
 Bioreactor vessels 50ml, 20pcs and sterile TPP TubeSpin®
 Bioreactor vessels, 50ml, with pH and O2 sensors, 10pcs
 Personal multi-channel bioreactor with non-invasive real time OD pH and pO2 measurement

RTS-8 plus is a personal bioreactor that utilizes patented Reverse-Spin® technology that applies non-invasive, mechanically driven, low energy consumption, innovative type of agitation where cell suspension is mixed by the ...

[read more](#)