


## External pre-filter set (polyphosphate/carbon/1 µm) with manometer

### DESCRIPTION



### CAT. NUMBER

**BS-070104-LK**

External pre-filter set (polyphosphate/carbon/1 µm) with manometer

### ACCESSORIES



**Labaqua Trace**  
BS-070105-A02  
ultrapure water system

Labaqua ultrapure systems are multi-purpose water purification systems. The Labaqua systems produce ultrapure and pure water directly from tap water.

Any configuration of a Labaqua ultrapure system produces both ultrapure ...

[read more](#)



**Labaqua HPLC**  
BS-070104-A02  
ultrapure water system

Labaqua ultrapure systems are multi-purpose water purification systems. The Labaqua systems produce ultrapure and pure water directly from tap water.

**Labaqua HPLC** produces water with very low organic carbon (TOC) ...

[read more](#)



**Labaqua Bio**  
30l storage tank  
ultrapure water system

Labaqua ultrapure systems are multi-purpose water purification systems. The Labaqua systems produce ultrapure and pure water directly from tap water.

**Labaqua Bio** system produces water with very low organic and ...

[read more](#)