


# Deionization module

## ■ DESCRIPTION



## ■ CAT. NUMBER

BS-070104-IK

Deionization module

## ■ ACCESSORIES



Labaqua Trace  
BS-070105-A02  
ultrapure water system



Labaqua HPLC  
BS-070104-A02  
ultrapure water system



Labaqua Bio  
30L storage tank  
ultrapure water system

Labaqua ultrapure systems are multi-purpose water purification systems. The Labaqua systems produce ultrapure and pure water directly from tap water.

Any configuration of a Labaqua ultrapure system produces both ultrapure ...

[read more](#)

[read more](#)

[read more](#)

Labaqua ultrapure systems are multi-purpose water purification systems. The Labaqua systems produce ultrapure and pure water directly from tap water.

Labaqua HPLC produces water with very low organic carbon (TOC) ...

Labaqua ultrapure systems are multi-purpose water purification systems. The Labaqua systems produce ultrapure and pure water directly from tap water.

Labaqua Bio system produces water with very low organic and ...