


# UV bulb 185 nm

## DESCRIPTION



## CAT. NUMBER

BS-070104-DK

UV bulb 185 nm

## ACCESSORIES



**Labagua Trace**  
BS-070105-A02  
ultrapure water system

Labagua ultrapure systems are multi-purpose water purification systems. The Labagua systems produce ultrapure and pure water directly from tap water.

Any configuration of a Labagua ultrapure system produces both ultrapure ...

[read more](#)



**Labagua HPLC**  
BS-070104-A02  
ultrapure water system

Labagua ultrapure systems are multi-purpose water purification systems. The Labagua systems produce ultrapure and pure water directly from tap water.

**Labagua HPLC** produces water with very low organic carbon (TOC) ...

[read more](#)



**Labagua Bio**  
30l storage tank  
ultrapure water system

Labagua ultrapure systems are multi-purpose water purification systems. The Labagua systems produce ultrapure and pure water directly from tap water.

**Labagua Bio** system produces water with very low organic and ...

[read more](#)