


Stacking kit for 3 × ES-20/80 | ES-20/80C

DESCRIPTION



CAT. NUMBER

BS-010167-PK

Stacking kit for 3× ES-20/80 | ES-20/80C

ACCESSORIES



ES-20/80

Software included, without platform
Orbital Shaker-Incubator

ES-20/80 shaker-incubator for biotechnological and pharmaceutical laboratories is a professional category equipment. The typical applications include - microbial and cell culture cultivation, protein expression, solubility studies, general mixing, as well ...

[read more](#)



ES-20/80C

Software included, without platform
Orbital Shaker-Incubator with cooling

ES-20/80C shaker-incubator for biotechnological and pharmaceutical laboratories is a professional category equipment. The typical applications include - microbial and cell culture cultivation, protein expression, solubility studies, general mixing, as well ...

[read more](#)