

BN-11/30, bucket set for rotor R-6P

DESCRIPTION

Spare bucket set for rotor R-6P. Plastic buckets (POM Kocetal) (6 pcs).

SPECIFICATIONS

Tube holders material	POM Kocetal
Tube's volume	50 ml
Max. speed	4200 rpm
Max. temperature	150 °C
Tube manufacturers	Falcon, Greiner Bio-one, Sarstedt, Corning, Nunc, etc.
Autoclavable	+



CAT. NUMBER

BS-010208-ZK

BN-11/30

ACCESSORIES



LMC-3000
 Without rotor
 Laboratory Centrifuge

LMC-3000 is a modern low-speed bench-top centrifuge designed for operation with microtest plates and centrifuge tubes up to 50 ml. This device is widely used in biomedical profile laboratories.

Features: ...

[read more](#)



LMC-12
 Without rotor
 Laboratory Centrifuge

LMC-3000 is a modern low-speed bench-top centrifuge designed for operation with microtest plates and centrifuge tubes up to 50 ml. This device is widely used in biomedical profile laboratories.

Features: ...

[read more](#)



LMC-4200R
 Without rotor
 Refrigerated Laboratory
 Centrifuge

Laboratory bench-top centrifuge with refrigeration **LMC-4200R** provides temperature control of biomaterial during centrifugation. Temperature control of the so-called "cold-shelf" is a gold standard for enzymologists and cell biologists because it ...

[read more](#)



R-6P
 BS-010208-XK
 Rotor

Swing-out rotor for conical bottom centrifuge tube. Plastic buckets (POM Kocetal).



Ri-6P
 BS-010221-IK
 Rotor

Swing-out rotor for conical bottom centrifuge tube. Plastic buckets (POM Kocetal).