

BN-16/90, bucket set for rotor R-12/10 R-24/10

DESCRIPTION

Spare bucket set for rotor R-12/10, R-24/10. Plastic buckets (12 pcs).

SPECIFICATIONS

Tube's volume	10–15 ml
Max. speed	4200 rpm
Tube manufacturers	Falcon, Greiner Bio-one, Sarstedt, Corning, Nunc, etc.
Autoclavable	+



CAT. NUMBER

BS-010208-SK	BN-16/90
--------------	----------



LMC-3000
Without rotor
Laboratory Centrifuge

LMC-3000 is a modern low-speed bench-top centrifuge designed for operation with microtest plates and centrifuge tubes up to 50 ml. This device is widely used in biomedical profile laboratories.

Features: ...

[read more](#)



LMC-12
Without rotor
Laboratory Centrifuge

LMC-3000 is a modern low-speed bench-top centrifuge designed for operation with microtest plates and centrifuge tubes up to 50 ml. This device is widely used in biomedical profile laboratories.

Features: ...

[read more](#)



LMC-4200R
Without rotor
Refrigerated Laboratory Centrifuge

Laboratory bench-top centrifuge with refrigeration **LMC-4200R** provides temperature control of biomaterial during centrifugation. Temperature control of the so-called "cold-shelf" is a gold standard for enzymologists and cell biologists because it ...

[read more](#)



R-12/10
BS-010208-BK
Rotor

Swing-out rotor for round bottom centrifuge tube.



Ri-12/10
BS-010221-MK
Rotor

Swing-out rotor for round bottom centrifuge tube.



R-24/10
BS-010212-JK
Rotor

Swing-out rotor for round bottom centrifuge tube.



Ri-24/10
BS-010221-LK
Rotor

Swing-out rotor for round bottom centrifuge tube.