

# Remote grade 2 water dispenser

## DESCRIPTION

Universal remote dispenser set with 3 m supply hose and water distribution module



## CAT. NUMBER

BS-070104-JK

Remote dispenser

## ACCESSORIES



**Labaqua Trace**  
BS-070105-A02  
ultrapure water system

Labaqua ultrapure systems are multi-purpose water purification systems. The Labaqua systems produce ultrapure and pure water directly from tap water.

Any configuration of a Labaqua ultrapure system produces both ultrapure ...

[read more](#)



**Labaqua HPLC**  
BS-070104-A02  
ultrapure water system

Labaqua ultrapure systems are multi-purpose water purification systems. The Labaqua systems produce ultrapure and pure water directly from tap water.

**Labaqua HPLC** produces water with very low organic carbon (TOC) ...

[read more](#)



**Labaqua Bio**  
30l storage tank  
ultrapure water system

Labaqua ultrapure systems are multi-purpose water purification systems. The Labaqua systems produce ultrapure and pure water directly from tap water.

**Labaqua Bio** system produces water with very low organic and ...

[read more](#)