

# TubeSpin® Bioreactor 50ml with pH and O2 sensors

## DESCRIPTION

Sterile TPP TubeSpin® Bioreactor vessel, 50 ml, with pH and O2 sensor spots, 1 pce.

## SPECIFICATIONS

Number of tubes in set	1
Volume	50 ml



## CAT. NUMBER

**BS-010170-AK**

Sterile TPP TubeSpin® Bioreactor vessel, 50 ml,  
with pH and O2 sensor spots, 1 pce.



#### RTS-8 Plus

Including TPP TubeSpin® Bioreactor vessels 50ml, 20pcs and sterile TPP TubeSpin® Bioreactor vessels, 50ml, with pH and O2 sensors, 10pcs  
Personal multi-channel bioreactor with non-invasive real time OD pH and pO2 measurement

RTS-8 plus is a personal bioreactor that utilizes patented Reverse-Spin® technology that applies non-invasive, mechanically driven, low energy consumption, innovative type of agitation where cell suspension is mixed by the ...

[read more](#)